

## TTTR

TTTR stands for Time-Tagged Time-Resolved. The basic concept of this measurement mode is to record every individual photon of a [TCSPC](#) experiment. Both the arrival time as well as the TCSPC timing information (i.e. the time from the preceding excitation pulse or to the following excitation pulse in reverse start stop experiments) are recorded.

The [TimeHarp 200](#) (also the [TimeHarp 100](#)), the [PicoHarp 300](#) [HydraHarp 400](#) support TTTR mode measurements. With the the [PicoHarp 300](#) a slight shift in the nomenclature was introduced. Due to its capability of simultaneously recording both signal inputs (in interactive mode used as 'Start' and 'Stop') as two separate TTTR traces with a time resolution of 4ps, a new TTTR mode was introduced, the so called T2-mode. The traditional recording as with the [TimeHarp 200](#) / [TimeHarp 100](#) is called [T3-mode](#). In [T2-mode](#) both signal inputs are functionally identical.

Today recording of [T3-mode](#) and [T2-mode](#) files is possible with all recent PicoQuant TCSPC devices ([TimeHarp 260](#), [PicoHarp 300](#) and [HydraHarp 400](#) and their included software, analysis is supported mainly by the [SymPhoTime](#) software package.

Due to its capability of including external TTL pulses as markers in the stream of photon events TTTR mode is suitable for [FLIM](#) measurements. The markers are used to synchronise the photon events with a scanning process allowing to sort the photons into image pixels.

Copyright of this document belongs to PicoQuant GmbH. No parts of it may be reproduced, translated or transferred to third parties without written permission of PicoQuant GmbH. All information given here is reliable to our best knowledge. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications and external appearances are subject to change without notice.



PicoQuant GmbH  
Rudower Chaussee 29 (IGZ)  
12489 Berlin  
Germany

P +49-(0)30-1208820-89  
F +49-(0)30-1208820-90  
info@picoquant.com  
www.picoquant.com